


# ADJUSTABLE SAFETY SHIELD WITH VELCRO STRAP AND THE FORHEAD FOAM PAD

## USER INSTRUCTION



The shield is made of transparent PVC 0,3-0,5 mm thick.  
Transparent PVC material is approved for food contact.

Shield flat dimensions is 300x400 mm.

The shield is resistant to alcohol-base disinfectants and can be reused after disinfection.

The shield doesn't adhere to the face (space approx. 3-4 cm). This distance allows the user to wear the hygienic mask underneath without risk of the fogging.

The soft foam pad protects the forehead against chafing and make the shield stable when using.

The velcro strap allows a comfortable fit of the mask to the head.

Shields should be stored at room temperature. Low temperature may cause cracks.

